






# GUIDELINE 21 CFR PART 11

How to implement the MAS-100 Atmos<sup>®</sup> into your compliant workflow



English



	OPTION 1	OPTION 2	OPTION 3
<b>21 CFR Part 11</b>	 <b>Out-of-the-box non 21 CFR Part 11 compliant</b>	<b>Stand-alone 21 CFR Part 11 supporting instrument</b>	<b>Smart link to a 21 CFR Part 11 compliant LIMS/EM software</b>
<b>User login for sampling</b>	 No user login on the instrument	21 CFR Part 11 compliant user management on the instrument	Logged in to LIMS/ EM software, but no user login on the instrument
<b>Data handling</b>	 If needed: Manual export of the tamper-proof audit trail	Manual export of the tamper-proof audit trail	Error-free and direct transfer of data via sampling barcode directly to LIMS/EM software
<b>Traceability</b>	 Limited traceability via Audit Trail Export: <ul style="list-style-type: none"> <li>- No user data available</li> <li>- Instrument data (Sampling Archive/ Audit Trail)</li> <li>- No sampling environment data (no scanner)</li> </ul>	Full traceability via Audit Trail Export: <ul style="list-style-type: none"> <li>- User data</li> <li>- Instrument data (Sampling Archive/ Audit Trail)</li> <li>- Sampling environment (scanner connected to instrument)</li> </ul>	Full traceability via external LIMS/EM software: <ul style="list-style-type: none"> <li>- Instrument data (in sampling barcode)</li> <li>- Sampling environment (scanner connected to LIMS/EM software)</li> </ul>
<b>Benefits</b>	 <ul style="list-style-type: none"> <li>- Out-of-the-box-operation</li> <li>- Maximum flexibility</li> <li>- For customers who work non-compliant (e.g. for customers outside pharma)</li> </ul>	<ul style="list-style-type: none"> <li>- Dedicated 21 CFR Part 11 compliant user management for MAS-100 Atmos</li> <li>- For customers who want to work compliant but do not use a LIMS/EM software</li> </ul>	<ul style="list-style-type: none"> <li>- Secure data transmission to higher-level systems</li> <li>- No additional user login required prior to sampling on the instrument</li> <li>- For customers using LIMS/EM software with connected peripheral devices (barcode scanner or tablet).</li> </ul>

**NOTE:** MAS-100 Atmos allows individual configuration and thus optimal adaptation even to use cases between the options shown above.