

Service Information Letter 7_2011

1 Release Note C&C PC SW 1.51 for MAS-100 NT / MAS-100 NT Ex

Instruments affected: C&C PC SW 1.50
Change level: optional
 recommend by MBV
 mandatory
 for information only

2 Reason

In the French language setting the user “calibreur” is not able to log in. The version 1.51 fixes this error.

There are no other changes to this version. All other languages are not affected and can work on 1.50 without restrictions.

3 Compatibility

Apart from this bug fix the version 1.50 is identical to V1.51 and therefore compatible.

4 Workaround

If you are on the French language setting on V1.50 and can't upgrade to V 1.51 you can either log in as “service” (this will also give access to the calibration setting) or change the language setting to any other language and use the calibrator account with the same password.

Stäfa, June 22nd 2011

Beat Glauser



MBV AG
CTO