

MAS-100 SIRIUS®


Microbial Air Sampler

QUICK START GUIDE




The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

FIRST STEPS



WIFI (OPTIONAL)




HOME BUTTON



1 Sec		On
<1 Sec		Menu
6 Sec		Reset

CONFIGURATION



Disable touch screen

Standby mode

Power off




Start measurement

Next page

IP address

Change password



Check lid integrity

For air lock or VHP

Next page


Button status





AIR. NOTHING ELSE.

Lot No:



Go to: certificate.mbv.ch



FILTER CERTIFICATE

Lot number archive

certificate.mbv.ch



DOWNLOADS

Quick tips, manuals and more

www.mbv.ch/sirius/files



SERVICE

Global MAS-100 service centers

www.mbv.ch/servicecenters



www.sigmaaldrich.com

The vibrant M and Millipore are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
LitCode: VM_UG14716EN V1 Date: 07/2025

